


3M™ Littmann® Classic II S.E. Stethoscope



3M Medical
from the 3M Health Care family

3M

3M™ Littmann® Classic II S.E. Stethoscope



Featuring a two-sided chestpiece, this classic stethoscope offers high acoustic sensitivity for superior performance. A tunable diaphragm on one side augments the traditional bell function of the opposite side.

Features

- Traditional bell combined with a tunable diaphragm that conveniently alternates between low- and high-frequency sounds without the need to turn over the chestpiece
- Excellent acoustic seal and comfortable fit with patented 3M™ Littmann® Snap Tight Soft-Sealing eartips
- Comfortably angled, anatomically correct headset
- Patient-friendly, nonchill rim and diaphragm
- Single-lumen tubing in a variety of popular colors
- Three-year warranty

Ordering Information

| Catalog # | Length | Tube Color |
|-----------|-----------------|----------------|
| 2201 | 28 in. (71,1cm) | Black |
| 2203 | 28 in. (71,1cm) | Gray |
| 2205 | 28 in. (71,1cm) | Navy Blue |
| 2206 | 28 in. (71,1cm) | Caribbean Blue |
| 2208 | 28 in. (71,1cm) | Hunter Green |
| 2209 | 28 in. (71,1cm) | Purple |
| 2210 | 28 in. (71,1cm) | Raspberry |
| 2211 | 28 in. (71,1cm) | Burgundy |
| 2812 | 28 in. (71,1cm) | Orange |
| 2813 | 28 in. (71,1cm) | Ceil Blue |
| 2814 | 28 in. (71,1cm) | Seafoam Green |

The Littmann stethoscopes listed do not contain natural rubber latex or dry natural rubber as components in the product or packaging.

For more information, call the 3M Health Care Helpline at **1-800-228-3957**, or visit us at **www.littmann.com**



Health Care

3M Center, Building 275-4W-02
St. Paul, MN 55144-1000
U.S.A.
1 800 228-3957

www.littmann.com

3M Canada Inc.

Post Office Box 5757
London, Ontario N6A4T1
Canada
1 800-563-2921

Outside of USA, please contact your 3M Representative



40% Pre-consumer waste paper
10% Post-consumer waste paper

Printed in U.S.A.
Copyright © 3M 2005
All Rights Reserved
70-2008-3506-7